

Product Name: Recombinant Human NMT1 Protein
Catalog #: PEH7538



Summary

Name	NMT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30419
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDCENEEDNS YNRGGLSPAN DTGAKKKKKK QKKKKEKGSE TDSAQDQPVK MNSLPAERIQ EIQAIELFS VGQGPAKTME EASKRSYQFW DTQPVPKLGE VVNTHGPVEP

image

SDS

Background

Alternative Names Glycylpeptide N-tetradecanoyltransferase 1 (EC 2.3.1.97) (Myristoyl-CoA:protein N-myristoyltransferase 1) (NMT 1) (Type I N-myristoyltransferase) (Peptide N-myristoyltransferase 1)

Note

For Research Use Only , Not for Diagnostic Use