


Product Name: Recombinant Human SRI Protein
Catalog #: PEH7540



Summary

Name	SORCN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P30626
Host	E.coli
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQDPLYGYFA AVAGQDGQID ADELQRCLTQ SGIAGGYKPF NLETCLMVS MLDRDMSGTM GFNEFKELWA VLNGWRQHFI SFDTDRSGTV DPQELQKALT TMGFRLSPQA VNSIAKRYST NGKITFDDYI ACCVKLRALT DSFRRRDTAQ

image



Background

Alternative Names Sorcin (22 kDa protein) (CP-22) (CP22) (V19)

Note

For Research Use Only , Not for Diagnostic Use