


Product Name: Recombinant Human DCTD Protein
Catalog #: PEH7545



Summary

Name	DCTD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P32321
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSEVSCKKRD DYLEWPEYFM AVAFLSAQRS KDPNSQVGAC IVNSENKIVG IGYNGMPNGC SDDVLPWRRT AENKLDTKYP YVCHAELENAI MNKNSTDVKG CSMYVALFPC NECAKLIQA GIKEVIFMSD KYHDSDEATA ARLLFNMAGV TFRKFIPKCS KIVIDFDSIN SRPSQKLQ

image



Background

Alternative Names Deoxycytidylate deaminase (EC 3.5.4.12) (dCMP deaminase)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DCTD Protein
Catalog #: PEH7545

