


Product Name: Recombinant Human USP6 Protein
Catalog #: PEH7546



Summary

Name	UBP6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35125
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QAQERKDILM KYDKGHRAGL PEDKGPEPVG INSSIDRFGI LHETELPPVT AREAKKIRRE MTRTSKWMEM LGEWETYKHS

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 6 (EC 3.4.19.12) (Deubiquitinating enzyme 6) (Proto-oncogene TRE-2) (Ubiquitin thioesterase 6) (Ubiquitin-specific-processing protease 6)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use