


Product Name: Recombinant Human ADD1 Protein
Catalog #: PEH7550



Summary

Name	ADDA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35611
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FASDGDSGTC SPLRHSFQKQ QREKTRWLNS GRGDEASEEG QNGSSPKSKT KWTKEDGHRT STSAVPNLFV PLNTNPKEVQ EMRNKIREQN LQDIKTAGPQ SQVLCGVVMD RSLVQGELVT ASKAIIEKEY QPHVIVSTTG PNPFTTLTDR ELEEYRREVE RKQKGSEENL DEAREQKEKS PPDQPAVPHP PPSTPIKLEE DLVPEPTTGD DSDAATFKPT LPDLSPDEPS EALGFPMLEK EEEAHRPPSP TEAPTEASPE PAPDPAPVAE EAAPSAVEEG AAADPGSDGS PGKSPSKKKK KFRTPSFLKK

image



Background

Alternative Names Alpha-adducin (Erythrocyte adducin subunit alpha)

Product Name: Recombinant Human ADD1 Protein
Catalog #: PEH7550



Note

For Research Use Only , Not for Diagnostic Use