

Product Name: Recombinant Human GRK3 Protein
Catalog #: PEH7551



Summary

Name	ARBK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35626
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VYLQKYPPPL IPPRGEVNAA DAFDIGSFDE EDTKGIKLLD CDQELYKNFP LVISERWQQE VTETVYEAVN ADTDKIEARK RAKNKQLGHE EDYALGKDCI MHGYMLKLG N PFLTQWQRRY FYLFPNRLEW RGEGESRQNL LTMEQILSVE ETQIKDKKCI LFRIKGGKQF VLQCESDPEF VQWKKELNET FKEAQRLLRR APKFLNKPRS GTVELPKPSL

image



Background

Alternative Names Beta-adrenergic receptor kinase 2 (Beta-ARK-2) (EC 2.7.11.15) (G-protein-coupled receptor kinase 3)

Note

Product Name: Recombinant Human GRK3 Protein
Catalog #: PEH7551



For Research Use Only , Not for Diagnostic Use