


Product Name: Recombinant Human CPOX Protein
Catalog #: PEH7552



Summary

Name	HEM6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P36551
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSLGRPEEEE DELAHRCSF MAPPVTDLGE LRRRPGDMKT KMELLILETQ AQVCQALAQV DGGANFSVDR WERKEGGGGI SCVLQDGCVF EKAGVSISVV HGNSLSEAAK QMRSRGKVLK TKDGKLPFCA MGVSIVHPK NPHAPTIHFN YRYFEVEEAD GNKQWWFGGG CDLTPTYLNQ EDVHFHRTL KEACDQHGPD LYPKFKKWCD DYFFIAHRGE RRGIGGIFD DLDSPSKEEV FRFVQSCARA VVPSYIPLVK KHCDDSFPTQ EKLWQQLRRG RYVEFNLLYD RGTKFGLFTP GSRIESILMS LPLTARWEYM HSPSENSKEA EILEVLRHPR DWVR

image



Background

Alternative Names Oxygen-dependent coproporphyrinogen-III oxidase, mitochondrial (COX)
(Coprogen oxidase) (Coproporphyrinogenase) (EC 1.3.3.3)

Product Name: Recombinant Human CPOX Protein
Catalog #: PEH7552



Note

For Research Use Only , Not for Diagnostic Use