


Product Name: Recombinant Human CSK Protein
Catalog #: PEH7554



Summary

Name	CSK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P41240
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSAIQAAWPS GTECIAKYNF HGTAEQDLPF CKGDVLTIVA VTKDPNWWYKA KNKVGREGII PANYVQKREG VKAGTKLSLM PWFHGKITRE QAERLLYPPE TGLFLVREST NYPGDYTLCV SCDGKVEHYR IMYHASKLSI DEEVYFENLM QLVEHYTSDA DGLCTRLIKP KVMEGTVAHQ DEFYRSGWAL NMKELKLLQT IGKGEFGDVM LGDYRGNKVA VKCIKNDATA QAFLAEASVM TQLRHSNLVQ LLGVIVEEKG GLYIVTEYMA KGSLVDYLRG RGRSVLGGDC LLKFSLDVCE AMEYLEGNF VHRDLAARNV LVSEDNVAKV SDFGLTKEAS STQDTGKLPV KWTAPEALRE KKFSTKSDVW SFGILLWEIY SFGRVYPRI PLKDVVPRVE KGYKMDAPDG CPPAVYEVK NCWHLDAAMR PSFLQLREQL EHIKTHELHL

image



Background

Product Name: Recombinant Human CSK Protein
Catalog #: PEH7554



Alternative Names

Tyrosine-protein kinase CSK (EC 2.7.10.2) (C-Src kinase) (Protein-tyrosine kinase CYL)

Note

For Research Use Only , Not for Diagnostic Use