


Product Name: Recombinant Human TEC Protein
Catalog #: PEH7557



Summary

| | |
|------------------------------|---|
| Name | TEC_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P42680 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 41 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | LIKRSQKKK TSPLNYKERL FVLTKSMLTY YEGRAEKKYR KGFIDVSKIK CVEIVKNDDG VIPCQNKYPF QVVDANTLY IFAPSPQSRD LWVKKLKEEI KNNNNIMIKY HPKFWTDGSY QCCRQTEKLA PGCEKYNLFE SSIRKALPPA PETKKRRPPP PIPEEEDNS EEIVVAMYDF QAAEGHDLRL ERGQEYLILE KNDVHWWRAR DKYGNEGYIP SNYVTGKKS NLDQYEWYCR NMNRSKAEQL LRSEDKEGGF MVRDSSQPGL YTVSLYTKFG GEGSSGFRHY HIKETTTPK KYLLAEKHAF GSIPEIIEYH KHNAAGLVTR |

image



Background

Alternative Names Tyrosine-protein kinase Tec (EC 2.7.10.2)

Product Name: Recombinant Human TEC Protein
Catalog #: PEH7557



Note

For Research Use Only , Not for Diagnostic Use