


Product Name: Recombinant Human LRPPRC Protein
Catalog #: PEH7559



Summary

Name	LPPRC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42704
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IAAKEKDIQE ESTFSSRKIS NQFDWALMRL DLSVRTGRI PKKLLQKVFN DTCRSGGLGG SHALLLRSC GSLLPELKE ERTEFAHRIW DTLQKLGAVY DVSHYNALLK VYLQNEYKFS PTDFLAKMEE ANIQPNRVTY QRLIASYCNV GDIEGASKIL GFMKTKDLPV TEAVFSALVT GHARAGDMEN AENILTVMRD AGIEPGPDY LALLNAYA EKVEKSELHL MDRDLLQIIF SFSKAGYPQY

image



Background

Alternative Names	Leucine-rich PPR motif-containing protein, mitochondrial (130 kDa leucine-rich protein) (LRP 130) (GP130)
--------------------------	---

Product Name: Recombinant Human LRPPRC Protein
Catalog #: PEH7559



Note

For Research Use Only , Not for Diagnostic Use