


**Product Name: Recombinant Human ACAA2 Protein**  
**Catalog #: PEH7560**



## Summary

<b>Name</b>	THIM_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P42765
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	11 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	KHKISREECD KYALQSQRW KAANDAGYFN DEMAPIEVKT KKGKQTMQVD EHARPQTTL QLQKLPPVFK

## image



## Background

<b>Alternative Names</b>	3-ketoacyl-CoA thiolase, mitochondrial (EC 2.3.1.16) (Acetyl-CoA acetyltransferase) (EC 2.3.1.9) (Acetyl-CoA acyltransferase) (Acyl-CoA hydrolase, mitochondrial) (EC 3.1.2.-) (EC 3.1.2.1) (EC 3.1.2.2) (Beta-ketothiolase) (Mitochondrial 3-oxoacyl-CoA thiolase) (T1)
--------------------------	--

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human ACAA2 Protein**  
**Catalog #: PEH7560**

---

