

## Summary

<b>Name</b>	UBP5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P45974
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	33 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SVLPTIRVPK AGDRVHKDEC AFSFDTPESE GGLYICMNTF LGFGKQYVER HFNKTGQRVY LHLRTRRPK EEDPATGTGD PPRKKPTRLA IGVEGGFDLS EEKFELDEDV KIVILPDYLE IARDGLGGLP DIVRDRVTSA VEALLSADSA SRKQEVQAWD GEVRQVSKHA FSLKQLDNPA RIPPCGWKCS KCDMRENLWL NLTDGSILCG RRYFDGSGGN NHAVEHYRET GYPLAVKLGT ITPDGADVYS YDEDDMVLDP SLAEHLSHFG IDMLKMQKTD

## image



## Background

<b>Alternative Names</b>	Ubiquitin carboxyl-terminal hydrolase 5 (EC 3.4.19.12) (Deubiquitinating enzyme 5) (Isopeptidase T) (Ubiquitin thioesterase 5) (Ubiquitin-specific-processing protease 5)
--------------------------	---

**Product Name: Recombinant Human USP5 Protein**  
**Catalog #: PEH7564**

---



**Note**

For Research Use Only , Not for Diagnostic Use