

Product Name: Recombinant Human NSF Protein
Catalog #: PEH7565



Summary

Name	NSF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P46459
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CPTDELSLTN CAVVNEKDFQ SGQHVIVRTS PNHRYTFTLK THPSVVPGSI AFSLPQRKWA GLSIGQEIEV SLYTFDKAKQ CIGTMTIEID FLQKKSIDSN PYDSDKMAAE FIQQFNNQAF SVGQQLVFSF NEKLFGLLVK DIEAMDPSIL KGEPATGKRQ KIEVGLVGN SQVAFEKAEN SSLNLIGKAK TKENRQSIIN PDWNFEKMGI GGLDKFSDI FRRAFASRVF PPEIVEQMGC

image



Background

Alternative Names Vesicle-fusing ATPase (EC 3.6.4.6) (N-ethylmaleimide-sensitive fusion protein) (NEM-sensitive fusion protein) (Vesicular-fusion protein NSF)

Note

Product Name: Recombinant Human NSF Protein
Catalog #: PEH7565



For Research Use Only , Not for Diagnostic Use