


Product Name: Recombinant Human ARCN1 Protein
Catalog #: PEH7568



Summary

Name	COPD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48444
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KITLTCGRDG GLQNMELHGM IMLRISDDKY GRIRLHVENE DKKGVQLQTH PNVDKKLFTA ESLIGLKNPE KSFPVNSDVG VLKWRLQTTE ESFIPLTINC WPSESGNGCD VNIEYELQED NLELNDVVIT IPLPSGVGAP VIGEIDGEYR HDSRRNTLEW CLPVIDAKNK SGSLEFSIAG QPNDFFPVQV SFVSKKNYCN IQVTKVTQVD GNSPVRFSTE TTFLVDKYEI

image



Background

Alternative Names Coatomer subunit delta (Archain) (Delta-coat protein) (Delta-COP)

Note

Product Name: Recombinant Human ARCN1 Protein
Catalog #: PEH7568



For Research Use Only , Not for Diagnostic Use