


Product Name: Recombinant Human MARCKSL1 Protein
Catalog #: PEH7572



Summary

Name	MRP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49006
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGSQSSKAPR GDVTAEEAAG ASPAKANGQE NGHVKSNGL SPKGEGESPP VNGTDEAAGA TGDAIEPAPP SQGAEAKGEV PPKETPKKKK KFSFKKPKL SGLSFKRNRK EGGDSSASS PTEEEQEQGE IGACSDEGTA QEGKAAATPE SQEPQAKGAE ASAASEEEAG PQATEPSTPS GPESGPTPAS AEQNE

image



Background

Alternative Names MARCKS-related protein (MARCKS-like protein 1) (Macrophage myristoylated alanine-rich C kinase substrate) (Mac-MARCKS) (MacMARCKS)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MARCKSL1 Protein
Catalog #: PEH7572

