

Product Name: Recombinant Human TUFM Protein
Catalog #: PEH7574



Summary

Name	EFTU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49411
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YVRDKPHVNV GTIGHVDHGK TTLTAAITKI LAEGGGAKFK KYEEIDNAPE ERARGITINA AHVEYSTAAR HYAHTDCPGH ADYVKNMITG TAPLDGCILV VAANDGPMPQ TREHLLARQ IGVEHVVVVYV NKADAVQDSE MVELVELEIR ELLTEFGYKG EETPVIVGSA LCALEGRDPE LGLKSVQKLL DAVDTYIPVP ARDLEKPFL PVEAVYSVPG RGTVVTGTLE RGILKKGDEC ELLGHSKNIR TVVTGIEMFH KSLRAEAGD NLGALVRGLK REDLRRGLVM VKPGSIKPHQ KVEAQVYILS KEEGGRHKPF VSHFMPVMFS LTWDMACRII LPPEKELAMP GEDLKFNLIL RQPMILEKGQ RFTLRDGNRT IGTGLVTNTL AMTEEEKNIK

image



Background

Alternative Names Elongation factor Tu, mitochondrial (EF-Tu) (P43)

Product Name: Recombinant Human TUFM Protein
Catalog #: PEH7574



Note

For Research Use Only , Not for Diagnostic Use