


Product Name: Recombinant Human RANBP2 Protein
Catalog #: PEH7577



Summary

Name	RBP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49792
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YIASVQGSTP SPRQKSMKGF YFAKLYYEAK EYDLAKKYIC TYINVQERDP KAHRFLGLLY ELEENTDKAV ECYRRSVELN PTQKDLVLKI AELLCKNDVT DGRAKYWLER AAKLFGSPA IYKLKEQLLD

image



Background

Alternative Names	E3 SUMO-protein ligase RanBP2 (EC 2.3.2.-) (358 kDa nucleoporin) (Nuclear pore complex protein Nup358) (Nucleoporin Nup358) (Ran-binding protein 2) (RanBP2) (p270)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RANBP2 Protein
Catalog #: PEH7577

