


Product Name: Recombinant Human SEP₁ Protein
Catalog #: PEH7578



Summary

Name	SPS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49903
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSTRESFNPE SYELDKSFRL TRFTELKGTG CKVPQDVLQK LLESLQENHF QEDEQFLGAV MPRLGIGMDT CVIPLRHGGL SLVQTTDIY PIVDDPYMMG RIACANVLS LYAMGVTECD NMLMLLGVS N KMTDRERDKV MPLIIQGFKD AAEEAGTSVT GGQTVLNPWI VLGAVATTVC QPNEFIMPDN AVPGDVLVLT KPLGTQVAVA VHQWLDIPEK WNKIKLVVTQ EDVELAYQEA MMNMARLNRT AAGLMHTFNA HAATDITGFG ILGHAQNLAK QQRNEVSFVI HNLPLVLAKMA AVSKACGNMF GLMHGTCPET SGLLICLPR EQAARFCAEI KSPKYGEGHQ AWIIGIVEKG NRTARIIDKP RIEVAPQVA TQNVNPTPGA TS

image



Background

Alternative Names Selenide, water dikinase 1 (EC 2.7.9.3) (Selenium donor protein 1)

Product Name: Recombinant Human SEPHS1 Protein
Catalog #: PEH7578



(Selenophosphate synthase 1)

Note

For Research Use Only , Not for Diagnostic Use