

**Product Name: Recombinant Human RENBP Protein**  
**Catalog #: PEH7584**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | RENBP_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | P51606   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 50 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.   |
| <b>Sequence</b>              | MEKERETLQA WKERVGQELD RVVAFWMEHS HDQEHGGFFT CLGREGRVYD<br>DLKYVWLQGR QVWMYCRLYR TFERFRHAQL LDKAAKAGGEF LLRYARVAPP<br>GKKCAFVLTR DGRPVKVQRT IFSECFYTM MNELWRATGE VRYQTEAVEM<br>MDQIVHWVQE DASGLGRPQL QGAPAAEPMA VPMMLLNLE QLGEADEELA<br>GKYAELGDC ARRILQHVQR DGQAVLENVS EGGKELPGCL GRQQNPGHTL<br>EAGWFLRHC IRKGDPELRA HVIDKFLLLP FHSGWDPDHG GLFYFQDADN<br>FCPTQLEWAM KLWWPHSEAM IAFLMGYSDS GDPVLLRLFY QVAEYFRQF<br>RDPEYGEWFG YLSREGKVAL SIKGGPFKGC FHVPRCLAMC EEMLGALLSR<br>PAPAPSPAPT PACRGAE |

## image



## Background

**Product Name: Recombinant Human RENBP Protein**  
**Catalog #: PEH7584**



---

**Alternative Names**

N-acylglucosamine 2-epimerase (AGE) (EC 5.1.3.8) (GlcNAc 2-epimerase) (N-acetyl-D-glucosamine 2-epimerase) (Renin-binding protein) (RnBP)

**Note**

For Research Use Only , Not for Diagnostic Use