


**Product Name: Recombinant Human RAP1GDS1 Protein**  
**Catalog #: PEH7588**



## Summary

<b>Name</b>	GDS1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P52306
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DKEDDITELK TGS DLMVLLL LGDESMQKLF EGGKGSVFQR VLSWIPSNNH QLQLAGALAI ANFARN DANC IHMVDNGIVE KLMDLLDRHV EDGNVTVQHA ALSALRN LAI PVINKAKMLS AGVTEAVLKF LKSEMPVQF KLLGTLRMLI DAQAEAAEQL GKNVKLVERL VEWCEAKDHA GVMGESNRL SALIRHSKSK DVIKTIVQSG GIKHLVTMAT SEHVIMQNEA LVALALIAAL ELGTAEKDLE SAKLVQILHR LLADERSAPE IKYNSMVLIC ALMGSECLHK EVQDLAFLDV VSKLRSHENK SVAQQASLTE

## image



## Background

**Alternative Names** Rap1 GTPase-GDP dissociation stimulator 1 (Exchange factor smgGDS) (SMG GDS protein) (SMG P21 stimulatory GDP/GTP exchange protein)

**Product Name: Recombinant Human RAP1GDS1 Protein**  
**Catalog #: PEH7588**

---



**Note**

For Research Use Only , Not for Diagnostic Use