


Product Name: Recombinant Human DAPK1 Protein
Catalog #: PEH7591



Summary

Name	DAPK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53355
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DYYDTGEELG SGQFAVVKKC REKSTGLQYA AKFIKKRRTK SSRRGVSRED IEREVSILKE IQHPNVITLH EYENKTDVI LILELVAGGE LFDFLAEKES LTEEEATEFL KQILNGVYYL HSLQIAHFDL KPENIMLLDR NVPKPRIKII DFGLAHKIDF GNEFKNIFGT PEFVAPEIVN YEPLGLEADM WSIGVITYIL LSGASPFLGD TKQETLANVS AVNYEFEDEY FSNTSALAKD FIRLLVKDP KKRMTIQDSL QHPWIKPKDT QQALSRKASA VNMEKFKKFA

image



Background

Alternative Names	Death-associated protein kinase 1 (DAP kinase 1) (EC 2.7.11.1)
--------------------------	--

Product Name: Recombinant Human DAPK1 Protein
Catalog #: PEH7591



Note

For Research Use Only , Not for Diagnostic Use