


Product Name: Recombinant Human ACLY Protein
Catalog #: PEH7592



Summary

Name	ACLY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53396
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKELLYKFC T TSAIQNRFK YARVTPDTDW ARLQDHPWL LSQNLVVKPD QLIKRRGKLG LVGVNLTLDG VKSWLKPRLG QEATVGKATG FLKNFLIEPF VPHSQAEFFY VCIYATREGD YVLFHHEGGV DVGDVDAKAQ KLLVGVDEKL NPEDIKKHLL VHAPEDKKEI LASFISGLFN FYEDLYFTYL EINPLVVTKD GVYVLDLAAK VDATADYICK VKWGDIEFPP PFGREAYPEE AYIADLDAKS

image



Background

Alternative Names ATP-citrate synthase (EC 2.3.3.8) (ATP-citrate (pro-S)-lyase) (ACL) (Citrate cleavage enzyme)

Note

Product Name: Recombinant Human ACLY Protein
Catalog #: PEH7592



For Research Use Only , Not for Diagnostic Use