


Product Name: Recombinant Human SMTN Protein
Catalog #: PEH7595



Summary

Name	SMTN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53814
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGALRKLLV TADLAERRRI RSAIRELQRQ ELEREEAALA SKRFAERQD NKENWLHSQQ REAEQRAALA RLAGQLES MN DVEELTALLR SAGEYEERKL IRAAIRRVRA QEIEAATLAG RLYSGRPNSG SREDSKGLAA HRLEQCEVPE REEQEQQAEV SKPTPTPEGT SQDVTTVTLL LRAPPGSTSS SPASPSSSPT PASPEPPLEP AEAQCLTAEV PGSPEPPSP PKTTSPEPQE SPTLPSTEGQ VVNKLKSGPK ETPAAQSPTR GPSDTKRADV AGPRPCQRS LSVLSPRQPAQ NRESTPLASG PSSFQRAGSV RDRVHKFTSD SPMAARLQDG TPQAALSPLT

image



Background

Alternative Names Smoothelin

Product Name: Recombinant Human SMTN Protein
Catalog #: PEH7595



Note

For Research Use Only , Not for Diagnostic Use