


Product Name: Recombinant Human ATXN1 Protein
Catalog #: PEH7597



Summary

Name	ATX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P54253
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAAGATTPSQ RSQLEAYSTL LANMGSLST PGHKAEEQQQ QQQQQQQHQ HQQQQQQQQQ QQQQQHLSRA PGLITPGSPP PAQQNQYVHI SSSPQNTGRT ASPPAIPVHL HPHQTMIPHT LTLGPPSQVV MQYADSGSHF VPREATKKAE SSRLQQAIQA KEVLNGEMEK SRRYGAPSSA DLGLGKAGGK SVPHPYESRH VVVHPSPSDY SSRDPSGVRA SVMVLPNSNT PAADLEVQQA THREASPSTL NDKSGHLHGK PGHRSYALSP HTVIQTTHSA

image



Background

Alternative Names Ataxin-1 (Spinocerebellar ataxia type 1 protein)

Product Name: Recombinant Human ATXN1 Protein
Catalog #: PEH7597



Note

For Research Use Only , Not for Diagnostic Use