


Product Name: Recombinant Human PMS1 Protein
Catalog #: PEH7598



Summary

Name	PMS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P54277
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKTAETDVLV NKVSSGKNY SNVDTSVIPF QNDMHNDESG KNTDDCLNHQ ISIGDFGYGH CSSEISNIDK NTKNAFQDIS MSNVSWENSQ TEYSKTCFIS SVKHTQSENG NKDHIDESGE NEEEAGLENS SEISADEWSR GNILKNSVGE NIEPVKILVP EKSLPCKVSN NNYPEQMN LNEDSCNKKS NVIDNKSGKV TAYDLLSNRV IKKPMASAL FVQDHRPQFL IENPKTSLED ATLQIEELWK TLSEEEKLKY EEKATKDLER

image



Background

Alternative Names PMS1 protein homolog 1 (DNA mismatch repair protein PMS1)

Product Name: Recombinant Human PMS1 Protein
Catalog #: PEH7598



Note

For Research Use Only , Not for Diagnostic Use