


Product Name: Recombinant Human VCP Protein
Catalog #: PEH7605



Summary

Name	TERA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55072
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLSTAILKQK NRPNRLIVDE AINEDNSVVS LSQPKMDELQ LFRGDTVLLK GKKRREAVCI VLSDDTCSDE KIRMNRVVRN NLRVRLGDVI SIQPCPDVKY GKRIHVLPID DTVEGITGNL FEVYLKPYFL EAYRPIRKGD IFLVRGGMRA VEFKVETDP SPYCIVAPDT VIHCEGEPIK REDEEESLNE

image



Background

Alternative Names Transitional endoplasmic reticulum ATPase (TER ATPase) (EC 3.6.4.6) (15S Mg(2+)-ATPase p97 subunit) (Valosin-containing protein) (VCP)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VCP Protein
Catalog #: PEH7605

