


**Product Name: Recombinant Human ADK Protein**  
**Catalog #: PEH7609**



## Summary

<b>Name</b>	ADK_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P55263
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	41 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MAAAEEPKP KKLKVEAPQA LRENILFGMG NPLLDISAVV DKDFLDKYSL KPNDQILAED KHKELFDELV KFKVEYHAG GSTQNSIKVA QWMIQQPHKA ATFFGCIGID KFGEILKRKA AEAHVDAHYY EQNEQPTGTC AACITGDNRS LIANLAAANC YKKEKHLDE KNWMLVEKAR VCYIAGFFLT VSPESVLKVA HHAENNRIF TLNLSAPFIS QFYKESLMKV MPYVDILFGN ETEAATFARE QGFETKDIKE IAKKTQALPK MNSKRQRIVI FTQGRDDTIM ATESEVTAF VLDQDQKEII DTNGAGDAFV GGFLSQLVSD KPLTECIRAG HYAASIIIRR TGCTFPEKPD FH

## image



## Background

**Alternative Names** Adenosine kinase (AK) (EC 2.7.1.20) (Adenosine 5'-phosphotransferase)

**Product Name: Recombinant Human ADK Protein**  
**Catalog #: PEH7609**

---



**Note**

For Research Use Only , Not for Diagnostic Use