

Product Name: Recombinant Human ADAR Protein
Catalog #: PEH7610



Summary

Name	DSRAD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55265
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQAWNQHSGV VRPDGHSQGA PNSDPSLEPE DRNSTSVSED LLEPFI AVSA QAWNQHSGVV RPDGHSQGSP NSDPGLEPED SNSTALEDP LEFLDMAEIK EKICDYLFNV SDSSALNLAK NIGLTKARDI NAVLIDMERQ GDVYRQGTP PIWHLTDKKR ERMQIKRNTN SVPETAPAAI PETKRNAEFL TCNIPTSNAS NNMVTTEKVE NGQEPVIKLE NRQEARPEPA

image

SDS

Background

Alternative Names	Double-stranded RNA-specific adenosine deaminase (DRADA) (EC 3.5.4.37) (136 kDa double-stranded RNA-binding protein) (p136) (Interferon-inducible protein 4) (IFI-4) (K88DSRBP)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use