


Product Name: Recombinant Human EIF3B Protein
Catalog #: PEH7612



Summary

Name	EIF3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55884
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLSQEIQIKQI KKDLKKYSKI FEQKDRLSQS KASKELVERR RTMMEDFRKY RKMAQELYME QKNERLELRG GVDTDDELDSN VDDWEEETIE

image



Background

Alternative Names	Eukaryotic translation initiation factor 3 subunit B (eIF3b) (Eukaryotic translation initiation factor 3 subunit 9) (Prt1 homolog) (hPrt1) (eIF-3-eta) (eIF3 p110) (eIF3 p116)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use