


Product Name: Recombinant Human AP1S2 Protein
Catalog #: PEH7613



Summary

Name	AP1S2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56377
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WRDLKIVYKR YASLYFCCAI EDQDNELITL EIIHRYVELL DKYFGSVCEL DIIFNFEKAY FILDEFLLGG EVQETSKKNV LKAIEQADLL QEEAETPRSV

image



Background

Alternative Names	AP-1 complex subunit sigma-2 (Adaptor protein complex AP-1 subunit sigma-1B) (Adaptor-related protein complex 1 subunit sigma-1B) (Clathrin assembly protein complex 1 sigma-1B small chain) (Golgi adaptor HA1/AP1 adaptin sigma-1B subunit) (Sigma 1B subunit of AP-1 clathrin) (Sigma-adaptin 1B) (Sigma1B-adaptin)
--------------------------	--

Note

Product Name: Recombinant Human AP1S2 Protein
Catalog #: PEH7613



For Research Use Only , Not for Diagnostic Use