


Product Name: Recombinant Human PER3 Protein
Catalog #: PEH7615




Summary

Name	PER3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56645
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDPPPSVTSQ RREEEKWEAQ SEGHPFITSR SSSPLQLNLL QEEMPRPSES PDQMRRNTCP QTEYCVTGNN GSESSPATTG ALSTGSPPRE NPSHTASAL STGSPPMKNP SHPTASALST GSPPMKNPSH PTASTLSMGL PPSRTPSHPT ATVLSTGSPP SESPSRTGSA ASGSSDSSIIY LTSSVYSSKI SQNGQQSQDV QKKETFPNVA EEPiWRMIRQ TPERILMTYQ VPERVKEVVL KEDLEKLESM RQQQPQFSHG QKEELAKVYN WIQSQTVTQE IDIQACVTCE NEDSADGAAT SCGQVLVEDS

image



SDS



WB

Background

Alternative Names Period circadian protein homolog 3 (hPER3) (Cell growth-inhibiting gene 13 protein) (Circadian clock protein PERIOD 3)

Product Name: Recombinant Human PER3 Protein
Catalog #: PEH7615



Note

For Research Use Only , Not for Diagnostic Use