



Summary

Name	RIPK4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57078
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEDLCEKPPD EVKETAHDLD VKSPPEPRSE VVPARLKRAS APTFDNDYSL SELLSQLDSG VSQAVEGPEE LSRSSSESKL PSSGSGKRLS GVSSVDSAFA SRGSLSLFE REPSTSDLGT TDVQKKKLVLD AIVSGDTSKL MKILQPQDVD LALDSGASLL HLAVEAGQEE CAKWLLLNNA NPNLSNRRGS TPLHMAVERR VRGVVELLLA RKISVNAKDE DQWTALHFAA QNGDESSTRLLLEKNASVNE VDFEGRTPMH VACQHGQENI VRILLRRGVD VSLQGKDAWL PLHYAAWQGH LPIVKLLAKQ PGVSVNAQTL DGRTPHLHAA QRGHYRVARI LIDLCSDVNV CSLLAQTPLH VAAETGHTST ARLLLHRGAG KEAMTSDGYT ALHLAARNGH LATVKLLVEE KADVLARGPL NQTALHLAAA HGHSEVVEEL VSADVIDLFD EQGLSALHLA AQGRHAQTVE TLLRHGAHIN LQSLKFQGGH GPAATLLRRS

image



Background

Product Name: Recombinant Human RIPK4 Protein
Catalog #: PEH7617



Alternative Names

Receptor-interacting serine/threonine-protein kinase 4 (EC 2.7.11.1)
(Ankyrin repeat domain-containing protein 3) (PKC-delta-interacting protein
kinase)

Note

For Research Use Only , Not for Diagnostic Use