


Product Name: Recombinant Human EPPK1 Protein
Catalog #: PEH7618



Summary

Name	EPIPL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P58107
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGFFDPNTHE NLTYVQLLR CVPDPDTGLY MLQLAGRGSA VHQLSEELRC ALRDARVTPG SGALQGQSVS VWELLYREV SEDRRQDLLS RYRAGTLTVE ELGATLTSL AQAAQARAE AEAGSPRPDP REALRAATME VKVGRLRGRA VPVWDVLASG YVSRAAREEL LAEFGSGTLD LPALTRRLTA IIEEAEAPG ARPQLQDAWR GPREGPAGR GDGDSGRSQR EGQGEGETQE AAAATAAARR QEQTLRDATM EVQRGQFQGR PVSVDVLFV SYLSEARRDE LLAQHAAGAL GLPDLVAVLT RVIEETEERL SKVSFRGLRR QVSASELHTS GILGPETLRD LAQGKTLQE

image



Background

Alternative Names	Epiplakin (450 kDa epidermal antigen)
--------------------------	---------------------------------------

Product Name: Recombinant Human EPPK1 Protein
Catalog #: PEH7618



Note

For Research Use Only , Not for Diagnostic Use