

**Product Name: Recombinant Human COPS2 Protein**  
**Catalog #: PEH7623**



## Summary

|                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | CSN2_HUMAN                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Accession #</b>           | P61201                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Species</b>               | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Molecular Mass</b>        | 52 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>              | MSDMEDDFMC DDEEDYDLEY SEDSNSEPNV DLENQYNSK ALKEDDPKAA<br>LSSFQKVLEL EGEKGEWGFK ALQMIKINF KLTNFPPEMMN RYKQLLTYYR<br>SAVTRNYSEK SINSILDYIS TSKQMDLLQE FYETTLEALK DAKNDRLWFK<br>TNTKLGLKLYL EREEYGKLQK ILRQLHQSCQ TDDGEDDLKK GTQLLEIYAL<br>EIQMYTAQKN NKKLKALYEQ SLHIKSAIPH PLIMGVIREC GGKMHLREGE<br>FEKAHTDFFE AFKNYDESGS PRRTTCLKYL VLANMLMKSG INPFDSQEAK<br>PYKNDPEILA MTNLVSAYQN NDITEFEKIL KTNHNSNIMDD PFIREHIEEL<br>LRNIRTQVLI KLIKPYTRIH IPFISKELNI DVADVESLLV QCILDNTIHG RIDQVNQLLE<br>LDHQKRGGAR YTALDKWTNQ LNSLNQAVVS KLA |

## image



## Background

**Product Name: Recombinant Human COPS2 Protein**  
**Catalog #: PEH7623**



---

**Alternative Names**

COP9 signalosome complex subunit 2 (SGN2) (Signalosome subunit 2) (Alien homolog) (JAB1-containing signalosome subunit 2) (Thyroid receptor-interacting protein 15) (TR-interacting protein 15) (TRIP-15)

**Note**

For Research Use Only , Not for Diagnostic Use