


**Product Name: Recombinant Human PSMC5 Protein**  
**Catalog #: PEH7626**



## Summary

<b>Name</b>	PRS8_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P62195
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	46 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LEEGKAGSGL RQYYLSKIEE LQLIVNDKSQ NLRRLQAQRN ELNAKVRLLR EELQLLQEQG SYVGEVVRAM DKKKVLVKVH PEGKFVVDVD KNIDINDVTP NCRVALRNDS YTLHKILPNK VDPLVSLMMV EKVPDSTYEM IGGLDKQIKE IKEVIELPVK HPELFEALGI AQPKGVLLYG PPGTGKTLA RAVAHHTDCT FIRVSGSELV QKFIGEGARM VRELFVMARE HAPSIIFMDE IDSIGSSRLE GGSGGDSEVQ RTMLELLNQL DGFEATKNIK VIMATNRIDI LDSALLRPGR IDRKIEFPPP NEEARLDILK IHSRKMNLTR GINLRKIAEL MPGASGAEVK GVCTEAGMYA LRERRVHVTQ EDFEMAVAKV MQKDSEKNMS

## image



## Background

**Alternative Names** 26S proteasome regulatory subunit 8 (26S proteasome AAA-ATPase subunit)

**Product Name: Recombinant Human PSMC5 Protein**  
**Catalog #: PEH7626**



---

RPT6) (Proteasome 26S subunit ATPase 5) (Proteasome subunit p45)  
(Thyroid hormone receptor-interacting protein 1) (TRIP1) (p45/SUG)

**Note**

For Research Use Only , Not for Diagnostic Use