


Summary

Name	UB2G1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P62253
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTELQSALLL RRQLAELNKN PVEGFSAGLI DDNDLYRWEV LIIGPPDTLY EGGVFKAHLT FPKDYPLRPP KMKFITEIWH PNVDKNGDVC ISILHEPGED KYGYEKPEER WLPIHTVETI MISVISMLAD PNGDSPANVD AAKEWREDRN GEFKRKVARC VRKSQETAFA

image



Background

Alternative Names	Ubiquitin-conjugating enzyme E2 G1 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme G1) (E217K) (UBC7) (Ubiquitin carrier protein G1) (Ubiquitin-protein ligase G1) [Cleaved into: Ubiquitin-conjugating enzyme E2 G1, N-terminally processed]
--------------------------	--

Product Name: Recombinant Human UBE2G1 Protein
Catalog #: PEH7627



Note

For Research Use Only , Not for Diagnostic Use