


Product Name: Recombinant Human RAB11A Protein
Catalog #: PEH7628



Summary

Name	RB11A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P62491
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGTRDDEYDY LFKVVLIGDS GVGKSNLLSR FTRNEFNLES KSTIGVEFAT RSIQVDGKTI KAQIWDTAGQ ERYRAITSAY YRGAVGALLV YDIAKHLTYE NVERWLKELR DHADSNIVIM LVGNKSDLRH LRAVPTDEAR AFAEKNGLSF IETSALDSTN VEAFAQTILT EIYRIVSQKQ MSDRRENDMS PSNNVPIHV PPTTENKPKV QCC

image



Background

Alternative Names Ras-related protein Rab-11A (Rab-11) (YL8)

Note

Product Name: Recombinant Human RAB11A Protein
Catalog #: PEH7628



For Research Use Only , Not for Diagnostic Use