


Product Name: Recombinant Human GNAI1 Protein
Catalog #: PEH7634



Summary

Name	GNAI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P63096
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGCTLSAEDK AAVERSKMID RNLREDGEKA AREVKLLLLG AGESGKSTIV KQMKIIHEAG YSEEECKQYK AVVYSNTIQS IIAIIRAMGR LKIDFGDSAR ADDARQLFVL AGAAEEGFMT AELAGVIKRL WKDSGVQACF NRSREYQLND SAAYYLNDLD RIAQPNYIPT QQDVLRTVK TTGIVETHFT FKDLHFKMFD VGGQRSEKWK WHCFEGVTA IIFCVASDY DLVLAEDEEM NRMHESMKLF DSICNNKWFT DTSILFLNK KDLFEKIKK SPLTICYPEY AGSNTYEEAA AYIQCCFEDL NKRKDTKEIY THFTCATDTK NVQFVFDAVT DVIKNNLKD CGLF

image



Background

Alternative Names	Guanine nucleotide-binding protein G(i) subunit alpha-1 (Adenylate cyclase-inhibiting G alpha protein)
--------------------------	--

Product Name: Recombinant Human GNAI1 Protein
Catalog #: PEH7634



Note

For Research Use Only , Not for Diagnostic Use