


Product Name: Recombinant Human PAFAH1B2 Protein
Catalog #: PEH7636



Summary

Name	PA1B2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P68402
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSQGDSNPAA IPHAAEDIQG DDRWMSQHNR FVLDCCKDEP DVLFVGDMSV QLMQQYEIWR ELFSPLHALN FGIGGDTTRH VLWRLKNGEL ENIKPKVIVV WVGTTNNHENT AEEVAGGIEA IVQLINTRQP QAKIIVLGLL PRGEKPNPLR QKNAKVNQLL KVSPLKLANV QLLDTDGGFV HSDGAISCHD MFDFLHLTGG GYAKICKPLH ELIMQLEET PEEKQTTIA

image



Background

Alternative Names	Platelet-activating factor acetylhydrolase IB subunit alpha2 (EC 3.1.1.47) (PAF acetylhydrolase 30 kDa subunit) (PAF-AH 30 kDa subunit) (PAF-AH subunit beta) (PAFAH subunit beta)
--------------------------	--

Product Name: Recombinant Human PFAH1B2 Protein
Catalog #: PEH7636



Note

For Research Use Only , Not for Diagnostic Use