


## Summary

<b>Name</b>	IGBP1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P78318
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	40 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MAAEDELQLP RLPETFETGR QLLDEVEVAT EPAGSRIVQE KVFKGLDLLE KAAEMLSQLD LFSRNEDLEE IASTDLKYL VPAFQGALTM KQVNPSKRDL HLQRAREHFI NYLTQCHCYH VAEFELPKTM NNSAENHTAN SSMAYPSLVA MASQRQAKIQ RYKQKKELEH RLSAMKSAVE SGQADDERVR EYLLHLQRW IDISLEEIES IDQEIKILRE RDSSREASTS NSSRQERPPV KPFILTRNMA QAKVFGAGYP SLPTMTVSDW YEQHRKYGAL PDQGIKAAP EEFKAAQQQ EEQEEKEED DEQTLHRARE WDDWKDTHPR GYGNRQNMG

## image



## Background

<b>Alternative Names</b>	Immunoglobulin-binding protein 1 (B-cell signal transduction molecule alpha 4) (Protein alpha-4) (CD79a-binding protein 1) (Protein phosphatase 2/4/6)
--------------------------	--

**Product Name: Recombinant Human IGBP1 Protein**  
**Catalog #: PEH7637**



---

regulatory subunit) (Renal carcinoma antigen NY-REN-16)

**Note**

For Research Use Only , Not for Diagnostic Use