


Product Name: Recombinant Human IFI35 Protein
Catalog #: PEH7640



Summary

Name	IN35_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P80217
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	<p>EVPKSLVSNL RIHCPLLAGS ALITFDDPKV AEQVLQQKEH TINMEECLRL VQVQPLELPM VTTIQMSSQL SGRRVLTGF PASLRLSEEE LLDKLEIFFG KTRNGGGDVD VRELLPGSVM LGFARDGVAQ RLCQIGQFTV PLGGQQVPLR VSPYVNGEIQ KAEIRSQPVP RSVLVLNIPD ILDGPELHDV LEIHFQKPTR GGGEVEALTV VPQQQQGLAV</p>

image



Background

Alternative Names Interferon-induced 35 kDa protein (IFP 35) (Ifi-35)

Note

Product Name: Recombinant Human IFI35 Protein
Catalog #: PEH7640



For Research Use Only , Not for Diagnostic Use