


Product Name: Recombinant Human HMGN5 Protein
Catalog #: PEH7641



Summary

Name	HMGN5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P82970
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DMRQEPKRRS ARLSAMLVPV TPEVKPKRTS SSRKMKT KSD MMEENIDTSA QAVAETKQEA VVEEDYNENA KNGEAKITEA PASEKEIVEV KEENIEDATE KGGEKKEAVA AEVKNEEDQ KEDEEDQNEE

image



Background

Alternative Names	High mobility group nucleosome-binding domain-containing protein 5 (Nucleosome-binding protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use