


Product Name: Recombinant Human PLCB2 Protein
Catalog #: PEH7642



Summary

Name	PLCB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q00722
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGPEGVDGRV RELKDRLELE LLRQGEEQYE CVLKRKEQHV AEQISKMMEL AREKQAAELK ALKETSENDT KEMKKKLETK RLRIQGMTK VTSDKMAQER LKREINNSHI QEVVQVIKQM TENLERHQEK LEEKQAACLE QIREMEKQFQ KEALAEYEAR MKGLEAEVKE SVRACLRTCF PSEAKDKPER ACECPPELCE QDPLIAKADA

image



Background

Alternative Names	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase beta-2 (EC 3.1.4.11) (Phosphoinositide phospholipase C-beta-2) (Phospholipase C-beta-2) (PLC-beta-2)
--------------------------	--

Product Name: Recombinant Human PLCB2 Protein
Catalog #: PEH7642



Note

For Research Use Only , Not for Diagnostic Use