



Product Name: Recombinant Human HNRNPU Protein
Catalog #: PEH7644



Summary

Name	HNRPU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q00839
Host	E.coli
Species	Human
Molecular Mass	75 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKVSELKEEL KKRRLSDKGL KAELMERLQA ALDDEEAGGR PAMEPGNGSL DLGGDSAGRS GAGLEQEAAA GGDEEEEEEE EEEEGISALD GDQMELGEEN GAAGAADSGP MEEEEAASED ENGDDQGFQE GEDELGDEEE GAGDENGHGE QQPQPPATQQ QQPQQQRGAA KEAAGKSSGP TSLFAVTVAP PGARQGQQQA GGKKKAEGGG GGGRPGAPAA GDGKTEQKGG DKRGVVRPR EDHGRGYFEY IEENKYSRAK SPQPPVEEED EHFDDTVVCL DTYNCDLHFK ISRDRLSASS LTMESFAFLW AGGRASYGVS KGKVCFEMKV TEKIPVRHLY TKDIDIHEVR IGWSLTTSKM LLGEEFSYG YSLKGIKTCN CETEDYGEKF DENDVITCFA NFESDEVELS YAKNGQDLGV AFKISKEVLA GRPLFPVLC HNCVEFNFG QKEKPYFPIP EETFIQNPV LEDRVRGPKG PEEKKDCEVV MMIGLPGAGK TTWVTKHAAE NPGKYNILGT NTIMDKMMVA GFKKQMA DTG KLNTLLQRAP QCLGKFIEIA ARKKRNFI LD QTNVSAAAQR RKMCLFAGFQ RKAVVVC PKD EDYKQRTQKK AEVEGKDLPE HAVLKMKG NF TLPEVAECFD EITYVELQKE EAQKLLEQYK

image



Product Name: Recombinant Human HNRNPU Protein
Catalog #: PEH7644



SDS

WB

Background

Alternative Names

Heterogeneous nuclear ribonucleoprotein U (hnRNP U) (GRIP120) (Nuclear p120 ribonucleoprotein) (Scaffold-attachment factor A) (SAF-A) (p120) (pp120)

Note

For Research Use Only , Not for Diagnostic Use