

**Product Name: Recombinant Human RELB Protein**  
**Catalog #: PEH7646**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | RELB_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | Q01201   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 37 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.   |
| <b>Sequence</b>              | EQPKQRGMRF RYECEGRSAG SILGESTEA SKTLPAIELR DCGGLREVEV<br>TACLVWKDWP HRVHPHSLVG KDCTDGICRV RLRPHVSPRH SFNNLGIQCV<br>RKKEIEAAIE RKIQLGIDPY NAGSLKNHQE VDMNVVRI CF QASYRDQQGQ<br>MRRMDPVLSE PVYDKKSTNT SELRICRINK ESGPCTGGEE LYLLCDKVQK<br>EDISVVSRA SWEGRADFSQ ADVHRQIAIV FKTPPYEDLE IVEPVTVNVF<br>LQRLTDGVCS EPLPFTYLPR DHDSYGVDKK RKRGMPPDVLG ELNSSDPHGI<br>ESKRRKKKPA |

## image



## Background

**Alternative Names** Transcription factor RelB (I-Rel)

**Product Name: Recombinant Human RELB Protein**  
**Catalog #: PEH7646**

---



**Note**

For Research Use Only , Not for Diagnostic Use