


Product Name: Recombinant Human PFKP Protein
Catalog #: PEH7650



Summary

Name	PFKAP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q01813
Host	E.coli
Species	Human
Molecular Mass	74 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CNVAVINVGA PAAGMNAAVR SAVRVGIADG HRMLAIYDGF DGFAGQIQKE IGWTDVGGWT GQGGSI LGTK RVLP GKYLEE IATQMRTHSI NALLIIGGFE AYLG LLELSA AREKHEEFCV PMVMVPATVS NNVP GSDFSI GADTALNTIT DTC DRIKQSA SGTKRRVFII ETMGGYCGYL ANMGGLAAGA DAA YIFEEP F DIRDLQSNVE HLTEKM KTTI QRGLVLRNES CSENYTTDFI YQLYSEEGKG VFDCRKNVLG HMQQGGAPSP FDRNFGTKIS

image



Background

Alternative Names	ATP-dependent 6-phosphofructokinase, platelet type (ATP-PFK) (PFK-P) (EC 2.7.1.11) (6-phosphofructokinase type C) (Phosphofructo-1-kinase isozyme C) (PFK-C) (Phosphohexokinase)
--------------------------	--

Product Name: Recombinant Human PFKP Protein
Catalog #: PEH7650



Note

For Research Use Only , Not for Diagnostic Use