


Product Name: Recombinant Human POU3F1 Protein
Catalog #: PEH7652



Summary

Name	PO3F1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03052
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DKIAAQGRKR KKRTSIEVGV KGALESHFLK CPKPSAHEIT GLADSLQLEK EVVRVWFCNR RQKEKRMTPA AGAGHPPMDD

image



Background

Alternative Names	POU domain, class 3, transcription factor 1 (Octamer-binding protein 6) (Oct-6) (Octamer-binding transcription factor 6) (OTF-6) (POU domain transcription factor SCIP)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use