



Product Name: Recombinant Human TNFAIP2 Protein
Catalog #: PEH7655



Summary

Name	TNAP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03169
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPLEAGAAPY REEEEAAKKK KEKKKSKGL ANVFCVFTKG KKKKGQPSSA EPEDAAGSRQ GLDGPPPTVE ELKAALERGQ LEAARPLLAL ERELAAAAAA GGVSEELVR RQSKVEALYE LLRDQVLGVL RRPLEAPPER LRQALAVVAE QEREDRQAAA AGPGTSGLAA TRPRRWLQLW RRGVAEAAEE RMGQRPAAGA

image



Background

Alternative Names Tumor necrosis factor alpha-induced protein 2 (TNF alpha-induced protein 2)
(Primary response gene B94 protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TNFAIP2 Protein
Catalog #: PEH7655

