



Product Name: Recombinant Human DNM1 Protein
Catalog #: PEH7656



Summary

Name	DYN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q05193
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IRKGWLTINN IGIMKGSKE YWFLTAENL SWYKDDEEKE KKYMLSVDNL KLRDVEKGFM SSKHIFALFN TEQRNVYKDY RQLELACETQ EEVDSWKASF LRAGVYPERV GDKEKASETE ENGSDSFMHS MDPQLERQVE TIRNLVDSYM AIVNKTVRDL MPKTIMHLM I NNTKEFIFSE LLANLYSCGD QNTLMEEAE QAQRRDEMLR MYHALKEALS IIGDINTTTV STPMPPPVDD SWLQVQSVPA GRRSPTSSPT PQRRAPAVPP ARPGSRGPAP GPPAGSALG GAPPVPSRPG ASDPFGPPP QVPSRPNRAP PGVPSRSGQA SPSRPESPRP

image



Background

Alternative Names Dynamin-1 (EC 3.6.5.5)

Product Name: Recombinant Human DNM1 Protein
Catalog #: PEH7656



Note

For Research Use Only , Not for Diagnostic Use