


Product Name: Recombinant Human MXD1 Protein
Catalog #: PEH7657



Summary

Name	MAD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q05195
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CDRKAVHQID QLQREQRHLK RQLEKLGIER IRMDSIGSTV SSERSDSDRE EIDVDVESTD YLTGDLDWSS SSVSDSDERG SMQSLGSDEG

image



Background

Alternative Names Max dimerization protein 1 (Max dimerizer 1) (Protein MAD)

Note

For Research Use Only , Not for Diagnostic Use